

DO NOT use this autoclave for the decontamination of biological and potentially infectious waste. Sterilization of glassware and liquids may be done in this particular autoclave.

- ♦ <u>Never put</u> sealed containers, solvents, volatile or corrosive chemicals, or radioactive materials in the autoclave.
- ◆ If the cycle did not meet the appropriate operating conditions or aborted, repeat the autoclave cycle. If the cycle fails again, report the malfunction to the responsible person listed below.
- ♦ Before unloading the Autoclave, <u>wait</u> for the pressure gauge to drop to zero with zero time remaining before opening the door.
- ◆ Using proper <u>heat resistant PPE</u>, carefully remove autoclaved materials after the cycle is completed and the appropriate waiting time for liquids has been observed.

Responsible Person for this Autoclave

Name:	
Phone:	



tp://ehs.uky.edu/ehs/biosafety ehsbiosafety@uky.edu